



Istem

MEDİKAL TIBBİ CİHAZ VE SANAYİ LTD.ŞTİ.

Intraosseous Infusion Needle



- Intraosseous Infusion Needle is designed to provide an intraosseous infusion when an intravenous route can not be established in an emergency, in patients with shock or loss of consciousness.
- Alternative way for fast, reliable and controlled access.
- Within ACLS (Advanced Cardiac Life Support), IV and central venous applications are kept at the same level.
- Drugs, blood and fluids arrive in seconds to venous system.
- Applicable by paramedic and first aid staff.
- Suitable for out-of-hospital use.
- With excellent tip design, the process is carried out in seconds without difficulty.

Product Name	Ref No	Diameter(G)	Length(CM)
Intraosseoz Infusion Needle-M Tip	ION1515	15	1,5
Intraosseoz Infusion Needle-M Tip	ION1525	15	2,5
Intraosseoz Infusion Needle-M Tip	ION1540	15	4,0
Intraosseoz Infusion Needle-M Tip	ION1545	15	4,5
Intraosseoz Infusion Needle-M Tip	ION1815	18	1,5



Istem

MEDİKAL TIBBİ CİHAZ VE SANAYİ LTD.ŞTİ.

Istem Medikal Tıbbi Cihaz ve San. Ltd. Şti.

Anadolu OSB. Mah. 29 Ekim Cad. No:39 Malıköy/Ankara/TURKEY

Phone: +90 312 394 55 62 Fax: +90 312 219 56 13

export@istemmedikal.com | www.istemmedikal.com